

Centrifuge

The worlds first 5ml microcentrifuge

This Mini Centrifuge is a benchtop centrifuge that is ideal for quick spin downs, microfiltration of HPLC samples and cell separations. It comes complete with a variety of rotors and features two operating speeds, three spin modes, and three choices of LED light colours. This is currently the only mini centrifuge on the market that accommodates 1.5ml and 5ml microtubes. With large centrifuge features packed into a personal benchtop unit, this is a cost-effective solution for laboratories on a budget.



LED colours to choose from



* Zisys Ltd reserve the right to change the design without prior notice as part of the ongoing design process.

Specifications

Speed	Slow: 6,000 RPM Fast: 10,000 RPM
Relative Centrifugal Force	Slow: 2,000 x G Fast: 4,800 x G
Noise Level	< 75dB
Rotor Capacity	Circular: 8 x 1.5ml tubes(8 x 0.5ml with supplied adapters) Strip: Two 8 x 0.2ml tube strips or 16 x 0.2ml tubes Microtube: 4 x 5ml microtubes
Speed Control	Variable
Power Requirements	100-240v ~ 1.0A; 50/60Hz (120v plug supplied) No Load: 6,000 RPM 16 Watts, 5.0A 10,000 RPM 18 Watts, 5.0A Fully Loaded: 6,000 RPM 26 Watts, 5.0A 10,000 RPM 40 Watts, 5.0A
Spin Modes	Continuous, Timed and Quick-Spin
Timer	1 to 15 minutes
Material	Lid: Clear polycarbonate Body & Rotors: Industrial ABS Resin Feet: Non-Skid natural rubber
Dimensions (LxWxH with lid closed)	9.45 x 7.36 x 5.47 in. (240 x 187 x 139mm)
Weight	3lb. 6oz. (1.6kg)
LED Colours	Red, Blue, Green

Exclusive 5ml Tubes



T2076 5ml Conical Microtube, Clear	250/Pk
T2077 5ml Conical Microtube, Blue	250/Pk
T2078 5ml Conical Microtube, Yellow	250/Pk
T2079 5ml Conical Microtube, Green	250/Pk
T2080 5ml Conical Microtube, Rose	250/Pk

5ml Microtubes also available in sterile packs of 25 10per box

T7246-901 1.5ml Conical Microtube, Clear	1000/PK
T7246-601 1.5ml Conical Microtube, Blue	1000/PK
T7246-401 1.5ml Conical Microtube, Yellow	1000/PK
T7246-501 1.5ml Conical Microtube, Green	1000/PK
T7246-201 1.5ml Conical Microtube, Rose	1000/PK

zisysTM

39 Mackintosh Place
South Newmoor Industrial Estate
Irvine
Ayrshire
Scotland KA11 4JT
Tel: +44 (0) 1294 219266
www.zisysautomation.com
allan@zisysautomation.com

* Zisys Ltd reserve the right to change the design without prior notice as part of the ongoing design process.