

Innovative, Flexible Tube Rack

More Than a Rack...

Introducing an innovative flexible pipetting work station. Its modular design allows researchers to create a rack based on their individual requirements. Starting with a tiered base stand; coloured inserts corresponding to individual tube sizes (0.5, 1.5, 5, 15 & 50mL) are snapped in place. Rack can accommodate 12 inserts. Inserts can be changed in seconds to accommodate changing sample requirements. Disposable pipetting basins can be placed in any row and will hold up to 25mL of reagent or solution.



They can be set side by side or "tiered" to allow for simultaneous viewing of all tubes. The small footprint (10 x 5 x 3 inches) and high tube capacity.



Changing inserts couldn't be easier. Simply pull the quick release tab and the insert lifts from the rack. Insert can be removed by grasping the finger tab located on the edge of each insert. To replace simply slip in the new insert, pull the quick release tab until insert snaps in place.



Each tube cavity of the colour coded insert is numbered to ease sample tracking and identification. Finger tabs are located on opposite ends to aid easy removal from the rack.



Inserts are also ideal for ice buckets. When placed on ice, the inserts create a stable and secure tube holder that can remain in the ice even when the tube is removed. Tubes stay dry and clean.



Inserts can be removed from the FlexiRack and used as a rack in a remote location for further processing.



The disposable polypropylene reagent basins will hold up to 25mL of reagent and features a large front edge to facilitate removal without spillage. The "V" shaped bottom ensures expensive reagents will not be left in the basin.

Features

Tiered Rack Design

Easy to view each tube in the rack.

Modular Design

Allows technician to create a rack that meets their individual requirements.

No need to have separate tube racks for various size tubes.

Coloured Tube Inserts

One rack can hold tubes from 0.5 - 50mL at the same time in any position on the rack.

Small Rack Footprint

Maximizes the number of tubes available without taking up valuable bench space.

Disposable Pipetting Basins

Convenient storage of solutions or reagents for pipetting into the rack

Polypropylene Construction

Inserts can be autoclaved up to +121^hC (including basins!)
Good chemical resistance.

Improved Productivity

FlexiRack can hold up to 144, 1.5mL tubes at one time!

Exclusive 5ml Tube! A Perfect Compliment to the FlexiRack....



Manufactured in polypropylene, this tube will withstand autoclaving (120°C for 15 minutes at 15psi), boiling and centrifugation to 14,000rpm. The 5ml tube can be spun in the same centrifuges that accommodates a 15ml tube. It is also suitable for freezing to -90°C. A frosted marking area on the tube cap is provided for sample tagging and identification. Available in both sterile or non-sterile.

The 5ml tube is the ideal alternative to 15ml conical tubes when using volumes less than 5ml. It's smaller size means you no longer need to insert your pipettor barrel inside the tube, only the tip needs to be inserted to draw out the contents.

R6000	Complete Rack System with brackets	1/EA
R6001	Rack without brackets	2/PK
R6002	Brackets	6/PK
R6003	PP 0.5ml Tube Insert, accommodates 12 tubes, red	5/PK
R6004	PP 1.5ml Tube Insert, accommodates 12 tubes, purple	5/PK
R6005	PP 5ml Tube Insert, accommodates 7 tubes, green	5/PK
R6006	PP 15ml Tube Insert, accommodates 4 tubes, blue	5/PK
R6007	PP 50ml Tube Insert, accommodates 3 tubes, orange	5/PK
R6008	PP V-Shaped Basin	50/PK
T2076	5ml Conical Microtube, clear	250/PK
T2077	5ml Conical Microtube, blue	250/PK
T2078	5ml Conical Microtube, yellow	250/PK
T2079	5ml Conical Microtube, green	250/PK
T2080	5ml Conical Microtube, rose	250/PK
T7246-901	1.5ml Conical Microtube, clear	1000/PK
T7246-601	1.5ml Conical Microtube, blue	1000/PK
T7246-401	1.5ml Conical Microtube, yellow	1000/PK
T7246-501	1.5ml Conical Microtube, green	1000/PK
T7246-201	1.5ml Conical Microtube, rose	1000/PK

5ml Microtubes also available in sterile packs of 25, 10per box



Ask for details regarding
The only centrifuge that
accomodates 0.5ml and 1.5ml tubes

zisysTM

39 Mackintosh Place
South Newmoor Industrial Estate
Irvine

Ayrshire
Scotland KA11 4JT

Tel: +44 (0) 1294 219266

www.zisys.co.uk

allan@zisys.co.uk